

Specific Surrogate Matrix Development Significantly Improves Insulin Measurements and Abrogates Cross-Reactivity With its Analogues

S. Fleire, M. Hyvönen, C. Köhler, F. Bernet, H. Faust, P. Brennecke, P. Struwe.
Celerion Switzerland AG, 8320 Fehraltorf, Switzerland



Introduction

Human insulin remains an important biopharmaceutical agent with new products regularly introduced in a steadily expanding market. Technical difficulties pose a challenge to the development of bioanalytical methods for measuring insulin (i.e., endogenous insulin, cross-reactivity with insulin analogues).


Here we demonstrate that the introduction of a specific surrogate matrix in an ELISA based bioassay for measuring human insulin, significantly suppressed the variability arising from endogenous insulin, typically observed in insulin depleted sera. Furthermore, a specially formulated dilution buffer included in the sample incubation step abrogates cross-reactivity with insulin analogues.

This novel insulin bioassay has been successfully validated where cross reactivity with a wide range of insulin analogues was completely abrogated. These results demonstrate the superiority of the novel bioassay, when compared with insulin depleted sera bioassays. Moreover, the assay has been successfully automated for future high through-put insulin measurements.

Assay Development

The developed insulin-specific assay used to quantify insulin in human serum samples is based on a commercial sandwich anti-human insulin colorimetric ELISA kit (Figure 1). The kit is based on two monoclonal antibodies that recognize distinct epitopes in human insulin: one for capture and a second for detection (HRP labeled).

Figure 1: Sandwich ELISA format




The assay is fast, simple, and robust, which allowed for the introduction of modifications aimed at circumventing technical difficulties commonly encountered during bioanalytical method development for insulin (i.e., endogenous insulin, cross-reactivity with insulin analogues).

Endogenous levels of insulin -ranging in normal individuals from 10 to 400 pM- represent a major challenge during bioanalytical assay development. Insulin-depleted sera employed as matrix for standard curve preparation show lot-to-lot variation in their response curves affecting long term studies.

To overcome this limitation, a buffer of defined composition was developed to replace insulin-depleted serum as the standard curve matrix (surrogate matrix buffer). To this end, dose-response curves of recombinant insulin prepared in different buffers were compared to a curve prepared in a validated insulin-depleted serum: the buffer with the resulting lowest mean deviation was selected as the surrogate matrix (Figure 2).

Figure 2: Comparison of calibration curves of human insulin prepared with two different insulin-depleted sera and with the surrogate matrix buffer.



An additional challenge in the development of insulin bioanalytical methods are the insulin analogues, as they may interfere with the accurate measurement of insulin. Initial tests showed that the assay cross-reacts particularly with two insulin analogues, Aspart and Glargine, while not with Degludec, Determir, Glulisine or Lispro.

Different strategies were explored to abrogate this interference, including using specific blocking anti-analogue antibodies, different washing buffers and sample dilution buffers. To test the different approaches, insulin-depleted serum samples were spiked with high concentrations of Aspart and Glargine insulin analogues. Best results were obtained with the inclusion of a dilution buffer during the sample incubation step (Table 1).

Table 1: Cross-reactivity of Aspart and Glargine insulin analogues with the insulin assay in the absence (W/O) or presence (W) of a sample dilution buffer.

| Insulin Analogue | % of cross-reactivity | |
|---------------------|-----------------------|-------------------|
| | W/O dilution buffer | W dilution buffer |
| Aspart (10000 pM) | 4.15 | < 0.15 |
| Glargine (10000 pM) | 16.74 | < 0.15 |

These two modifications allowed for the development of a highly specific bioanalytical insulin assay with an analytical range of 20 pM to 600 pM. Method validation confirmed the assay performance with the included modifications (see Assay Validation).

Assay Validation

The human insulin assay performance was validated following international standards. The standard curve was prepared in surrogate matrix buffer with an analytical range of 20 pM to 600 pM and quality control (QC) samples were prepared in low insulin level pool ([insulin] ≈ 20 pM). The following concentrations were selected: 20.9 pM (LLOQ-QC, corresponding to the un-spiked pool), 55 pM (LQC), 225 pM (MQC), 420 pM (HQC) and 600 pM (ULOQ-QC).

The following assay parameters were evaluated: precision and accuracy (standards and QC samples), parallelism, dilution linearity, selectivity, cross-reactivity with insulin analogues and stability. Additionally, the method precision and accuracy performance was validated on an automated system.

Precision and Accuracy

Precision and accuracy was evaluated in 27 (twenty-seven) experiments for standards, and in 6 (six) independent experiments (5 replicates of each control level / experiment) for quality controls.

Inter-run precision for standards

| | STD 1 (20 pM) | STD 2 (35 pM) | STD 3 (60 pM) | STD 4 (100 pM) | STD 5 (160 pM) | STD 6 (250 pM) | STD 7 (450 pM) | STD 8 (600 pM) |
|------|---------------|---------------|---------------|----------------|----------------|----------------|----------------|----------------|
| Mean | 19.5 | 35.1 | 60.1 | 98.7 | 163 | 248 | 450 | 600 |
| %CV | 7.3 | 3.6 | 2.4 | 2.2 | 2.2 | 2.0 | 1.4 | 0.6 |
| %Dev | -2.5 | 0.3 | 0.2 | -1.3 | 1.9 | -0.8 | 0.0 | 0.0 |
| n | 26 | 26 | 27 | 27 | 27 | 27 | 27 | 26 |

Inter-run precision for QC samples

| | LLOQ-QC (20.9 pM) | LQC (55.0 pM) | MQC (225 pM) | HQC (420 pM) | ULOQ-QC (600 pM) |
|------|-------------------|---------------|--------------|--------------|------------------|
| Mean | 21.2 | 55.2 | 222 | 421 | 594 |
| %CV | 5.5 | 4.9 | 4.1 | 4.4 | 5.0 |
| %Dev | 1.4 | 0.4 | -1.3 | 0.2 | -1.0 |
| n | 30 | 30 | 30 | 30 | 30 |

Parallelism

Parallelism was evaluated by diluting three individual samples with high insulin levels with surrogate matrix buffer (four dilutions within the analytical range were evaluated).

| Dilution Factor | Individual 1 | | | Individual 2 | | | Individual 3 | | |
|-----------------|----------------|-----|------|----------------|-----|------|----------------|-----|------|
| | [Insulin] (pM) | %CV | %Dev | [Insulin] (pM) | %CV | %Dev | [Insulin] (pM) | %CV | %Dev |
| Undiluted | 158 | 1.4 | | 172 | 0.6 | | 212 | 0.5 | |
| 2 | 182 | 0.5 | 15.2 | 188 | 0 | 9.3 | 218 | 1.1 | 2.8 |
| 3 | 189 | 2.5 | 19.6 | 192 | 0.9 | 11.6 | 224 | 0.2 | 5.7 |
| 4 | 186 | 0.7 | 17.7 | 201 | 0.5 | 16.9 | 225 | 0 | 6.1 |
| 6 | 188 | 0.8 | 19.0 | 203 | 3 | 18.0 | 220 | 1 | 3.8 |
| Mean | 180.6 | | | 191.2 | | | 219.8 | | |
| SD | 12.9 | | | 12.4 | | | 5.2 | | |
| %CV | 7.2 | | | 6.5 | | | 2.4 | | |

Dilution Linearity

Dilution linearity was evaluated with low insulin level pool sample spiked with 20000 pM of recombinant insulin. Serial dilutions were then performed with surrogate matrix buffer and the samples analysed.

| Dilution Factor | 0 | 10 | 40 | 100 | 400 |
|-----------------------------|-------|-------|------|------|------|
| Nominal Concentration (pM) | 20000 | 2000 | 500 | 200 | 50 |
| Measured concentration (pM) | >ULOQ | >ULOQ | 533 | 227 | 55.8 |
| SD | | | 11.7 | 5.12 | 1.50 |
| %CV | N/AP | N/AP | 2.2 | 2.3 | 2.7 |
| %Dev | | | 6.6 | 13.5 | 11.6 |
| n | 5 | 5 | 5 | 5 | 5 |

Selectivity

Twelve samples were evaluated for selectivity, including a lipemic sample. Selected samples were spiked with an additional 30 pM of recombinant insulin and analysed.

| Sample ID | Condition | Basal Level | Spike (+30.0 pM) | | |
|---|-----------|---------------------|------------------------|---------------------|------------|
| | | Measured Conc. (pM) | Theoretical Conc. (pM) | Measured Conc. (pM) | % Recovery |
| CM/17-0373 | | 39.8 | 69.8 | 64.9 | 93.0 |
| CM/17-0379 | | 55.2 | 85.2 | 80.8 | 94.8 |
| CM/17-0382 | | 41.8 | 71.8 | 69.1 | 96.2 |
| CM/17-0389 | | 26.9 | 56.9 | 57.0 | 100.2 |
| CM/17-0394 | | 31.9 | 61.9 | 62.0 | 100.2 |
| CM/17-0402 | | 30.3 | 60.3 | 55.1 | 91.4 |
| CM/17-0408 | | 30.3 | 60.3 | 58.9 | 97.7 |
| CM/17-0414 | | 20.7 | 50.7 | 49.8 | 98.2 |
| CM/17-0417 | | 29.1 | 59.1 | 58.7 | 99.3 |
| CM/17-0969 | | 24.8 | 54.8 | 51.5 | 94.0 |
| CM/17-0974 | | 40.2 | 70.2 | 66.3 | 94.4 |
| CM/17-1429 | lipemic | 86.9 | 117 | 110 | 94.1 |
| Number of serum samples assessed | | | | | 12 |
| Number of serum samples meeting acceptance criteria | | | | | 12 |
| % of serum samples meeting acceptance criteria | | | | | 100 |

Cross-Reactivity

Cross-reactivity of the assay with insulin analogues was evaluated during validation as the effect of the analogues on the measurement of both, the LLOQ-QC and ULOQ-QC.

Samples containing insulin at the LLOQ and ULOQ levels were spiked with different insulin analogues and the samples analysed. The deviation in insulin concentration was then evaluated.

The concentrations of analogues were selected based on previous PK studies analyzed at Celerion. They comprise in all cases the C_{max} values expected during therapeutics treatment with the corresponding analogues.

| | Aspart (2000 pM) | Degludec (10000 pM) | Determir (10000 pM) | Glargine (800 pM) | Glulisine (2000 pM) | Lispro (2000 pM) |
|---------------------|------------------|---------------------|---------------------|-------------------|---------------------|------------------|
| Nominal conc. (pM) | 20.9 | | | | | |
| Measured conc. (pM) | 22.6 | 22.4 | 21.2 | 21.8 | 21.5 | 19.7 |
| SD | 0.666 | 0.460 | 0.678 | 0.559 | 0.371 | 0.865 |
| %Dev | 8.1 | 7.2 | 1.4 | 4.3 | 2.9 | -5.7 |
| n | 5 | 5 | 5 | 5 | 5 | 5 |
| Nominal conc. (pM) | 600 | | | | | |
| Measured conc. (pM) | 573 | 553 | 576 | 567 | 559 | 559 |
| SD | 7.85 | 12.9 | 4.47 | 6.02 | 14.2 | 8.53 |
| %Dev | -4.5 | -7.8 | -4.0 | -5.5 | -6.8 | -6.8 |
| n | 5 | 5 | 5 | 5 | 5 | 5 |

Stability

Stability studies of biopharmaceuticals are routinely performed during method validation at Celerion. Short-term, long-term, bench-top and freeze/thaw cycles (-20°C and -80°C) of quality control samples are evaluated, as well as specific conditions such as stability in hemolytic samples. These studies are drug-specific, and therefore, they are not reported for this method.

Stability of standards in surrogate matrix buffer was also investigated in this case, as the newly developed buffer is a novel introduction for the assay. Standards were subjected to 4 freeze/thaw cycles at -20°C or -80°C and analyzed with a freshly prepared standard curve (only standards 1 and 8 are shown).

| Condition | STD 1 (20 pM) | | STD 8 (600 pM) | |
|-------------------|---------------|------|----------------|------|
| | Conc. (pM) | %Dev | Conc. (pM) | %Dev |
| Freeze/thaw: | 19.6 | -2.0 | 597 | -0.5 |
| 4 cycles at -20°C | 19.9 | -0.5 | 599 | -0.2 |
| | 20.8 | 4.0 | 593 | -1.2 |
| | 20.1 | 0.5 | 604 | 0.7 |
| | 21.1 | 5.5 | 597 | -0.5 |
| | 20.3 | 1.5 | 581 | -3.2 |
| Mean | 20.3 | | 595 | |
| %CV | 2.8 | | 1.3 | |
| %Dev | 1.5 | | -0.8 | |
| n | 6 | | 6 | |
| Freeze/thaw: | 21.1 | 5.5 | 625 | 4.2 |
| 4 cycles at -80°C | 20.7 | 3.5 | 587 | -2.2 |
| | 20.1 | 0.5 | 600 | 0.0 |
| | 20.6 | 3.0 | 604 | 0.7 |
| | 20.2 | 1.0 | 602 | 0.3 |
| | 20.2 | 1.0 | 604 | 0.7 |
| Mean | 20.5 | | 604 | |
| %CV | 1.9 | | 2.0 | |
| %Dev | 2.5 | | 0.7 | |
| n | 6 | | 6 | |

Automation

The assay precision and accuracy performance was also validated on an automated system. Each automated batch analysis comprises 8 plates; for validation purposes, plates 1 and 8 were used as test plates, while plates 2 to 7 were used as dummy plates (only inter-run precision and accuracy evaluation of quality control samples are shown).

| | LLOQ-QC (20.9 pM) | LQC (55.0 pM) | MQC (225 pM) | HQC (420 pM) | ULOQ-QC (600 pM) |
|------|-------------------|---------------|--------------|--------------|------------------|
| Mean | 19.2 | 54.7 | 206 | 362 | 495 |
| %CV | 9.0 | 3.9 | 4.3 | 5.9 | 3.8 |
| %Dev | -8.1 | -0.5 | -8.4 | -13.8 | -17.5 |
| n | 10 | 14 | 19 | 19 | 15 |

The automated method allows the analysis of ≈500 samples per day.

Conclusion

Bioanalytical assays for insulin are negatively impacted by the presence of endogenous insulin, as well as by the interference of insulin analogues. In the present work, a bioanalytical assay to measure the concentration of human insulin in human serum samples was developed and validated. The assay is based on a commercial kit and has an analytical range from 20 pM to 600 pM.

The two novel features introduced in this assay, i.e. use of surrogate matrix buffer and sample dilution buffer, significantly reduced the variability due to lot-to-lot variation of insulin-depleted sera and abrogated the cross-reactivity with insulin analogues.

Bioanalytical assays developed at Celerion Switzerland AG offer important tools for early clinical development of biopharmaceutical, as well as for the evaluation of biosimilars.