

Development and Optimization of a Highly Specific and Sensitive Assay to Detect Acetylcholinesterase Activity and Its Inhibition by Donepezil

Poster Number
04R0800

J. Thompson¹, R. Islam¹, M. Balasurya² and J. Knight²
¹Celerion and ²Somru BioScience



Background

Acetylcholinesterase (AChE) is one of the most critical enzymes for nerve response and function. It is primarily found at neuromuscular junctions and cholinergic synapses in the central nervous system. Cholinesterases are enzymes responsible for cleaving cholinergic neurotransmitters; this action allows cholinergic neurons to return to their resting state after activation. If cholinergic neurons were to remain activated, then muscles could not relax once they had been tensed. AChE activity serves to terminate the synaptic transmission. Vertebrates have two cholinesterases: acetylcholinesterase (AChE) is located mainly in chemical synapses and erythrocyte membranes while butyrylcholinesterase (BChE) is located mainly in plasma. AChE and BChE are differentiated based on their preferred substrate: AChE hydrolyzes acetylcholine (ATC) faster, while BChE hydrolyzes butyrylcholine (BTC) faster. However, both substrates can be hydrolyzed by both cholinesterases. This makes the specific measurement of AChE challenging.

Purpose and Methods

AChE inhibitors are among the key drugs approved for Alzheimer's disease (AD) and myasthenia gravis. In order to evaluate pharmacodynamic properties of these drugs it is critical to measure the specific activity of AChE without interference from BChE. Here we present an enzyme activity method based on the kit from Somru BioScience Inc. that is highly specific to AChE activity and sensitive enough to measure inhibition of AChE activity by donepezil.




Figure 1: Eel acetylcholinesterase (AChE) and equine butyrylcholinesterase (BChE) activities were tested with or without BChE inhibitor. a) Kinetic calibration curves show diminished OD values for BChE but not AChE, proving specificity of the inhibitor. b) Inhibition is concentration dependent: STD H kinetic curves are shown for AChE, in which activity remains uninhibited, and BChE, in which activity is inhibited at both 1:8 and 1:4 dilutions of inhibitor.




Figure 2: Contribution of background signal from whole blood and determination of optimal dilution for use in assay. Eel acetylcholinesterase (AChE) activity in human whole blood diluted 1:100, 1:200, or 1:400 was tested with or without substrate- dotted lines indicate background signal.

Results from 343 nm readout (Figures 1 and 2)

- Confirmed specificity of assay substrate is for AChE and not BChE
- High background noise, presumably from contribution of hemoglobin, necessitates high dilution of sample, possibly leading to loss of signal from AChE activity
- New chromophore using 670 nm wavelength will be tested




Figure 3: Area under curve for AChE activity in RBCs treated with donepezil at concentrations ranging from 7.5 – 150 ng/mL. Undiluted RBCs were analyzed at 670 nm in a kinetic assay with reads every 30 seconds. 2 minutes was chosen for subsequent method qualification and sample analysis.




Figure 4: Inhibition of AChE activity in RBCs by donepezil concentrations ranging from 7.5 – 150ng/mL. Undiluted RBCs were analyzed at 670 nm in a kinetic assay with reads every 30 seconds.




Figure 5: Sample stability in matrix. Inhibition of AChE activity in RBC was measured for various levels of donepezil treatment in samples assayed immediately after removal from -80°C, or in those that were allowed to sit at 50C for 24 hours. Assays were performed at 670 nm 2 minutes post-addition of substrate.




Figure 6: Inhibition of AChE activity in RBC correlates to concentration of donepezil. Samples from two individuals treated with donepezil were assayed at 670 nm for AChE activity 2 minutes post-addition of substrate. Assay is sensitive enough to detect inhibition of AChE activity at 7.5 ng/mL donepezil.

	untreated			7.5 ng/mL			25 ng/mL			50 ng/mL			100 ng/mL			150 ng/mL		
Well	mean	%CV	% inhibition	mean	%CV	% inhibition	mean	%CV	% inhibition	mean	%CV	% inhibition	mean	%CV	% inhibition	mean	%CV	% inhibition
81	85	81	83	71	69	70	70	66	68	49	47	48	52	53	53	50	50	50
82	85	79	82	76	73	75	65	65	65	46	47	46	51	52	51	48	47	47
83	85	82	84	76	70	73	65	67	66	48	45	47	49	48	48	44	44	47
mean	85	81	83	76	70	73	65	67	66	48	45	47	50	51	51	47	46	47
%CV	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3
background subtracted	59	49	42	42	42	42	42	42	42	23	23	23	23	23	23	23	23	23
% inhibition	18	18	18	18	18	18	18	18	18	61	61	61	61	61	61	61	61	61

Table 1: Data for 6 runs are shown in which 3 sets of untreated RBC samples and 3 sets of 5 levels of donepezil fortified RBC samples were assayed in duplicate. Acceptance criteria (% CV for all levels of donepezil treatment must be equal to or less than 20.0%) was met, showing assay precision.

Conclusion

- Use of substrate coupled with 670 nm chromophore abrogates background noise, so no need to dilute RBC samples
- High assay reproducibility
- Samples can be re-assayed after freeze thaw cycles
- Demonstration of short-term stability in matrix
- Inhibition of AChE activity directly correlates to donepezil concentration